



INTERNATIONAL DIAGNOSTIC & HEALTHCARE CENTER

Lab Admin Web

User Manual

A step-by-step guide for lab staff using the web administration dashboard: orders, sample collection, lab results, catalog, customers, pricing rules, reports and system settings.

Version 1.0 · September 2026

For authorized clinical and administrative staff only.

Contents

- 1. Introduction 3**
 - Who should read this manual 3
 - How to read this manual 3
- 2. Getting Started 4**
 - Signing in 4
 - The screen layout 4
 - Your profile and password 4
 - Language and theme 4
 - Notifications 5
 - Working with lists and tables 5
 - Signing out 5
- 3. The Order Lifecycle 6**
 - Priority 6
 - Report delivery method 6
 - Payment status 6
- 4. Orders 7**
 - The order list 7
 - Creating an order 7
 - Viewing order details 7
 - Adding tests to a pending order 8
 - Recording a payment 8
 - Deleting an order 8
- 5. Collection & Routing 9**
 - Planning a route 9
 - Assigning collectors 9
 - Updating a scheduled pickup 9
 - Starting collection 9
- 6. Lab Results 10**
 - Selecting an order 10
 - Step 1 — Receive the sample 10
 - Step 2 — Finish lab processing 10
 - Step 3 — Upload result PDFs 10
 - Step 4 — Run AI review 10
 - Step 5 — Release or deliver 10

7. Lab Tests (Test Catalog)	12
Adding a test	12
Editing, deactivating and deleting	12
8. Users	13
Adding a user	13
Editing a user	13
9. Staff	14
Adding a staff member	14
Editing, deactivating and resetting passwords	14
10. Ratings & Feedback	15
11. Discounts	15
Adding a discount	15
12. Referral Fees	16
Settings tab	16
Report tab	16
13. Membership Tiers	16
Adding a tier	16
14. Service Geofences	17
Adding a zone	17
Testing with the fee calculator	17
15. Advertisements	18
16. Loyalty Points	18
Adding an earn rule	18
Setting the redemption rate	18
17. System Settings	19
18. Reports & Analytics	19
19. Permissions (Admin only)	20
20. Troubleshooting & FAQ	20
Security reminders	20

1. Introduction

What the Lab Admin Web is, who uses it, and how this manual is organised.

The **Lab Admin Web** is the browser-based back office for the laboratory. Patients and partners place orders through the mobile app; lab staff use this website to review those orders, plan sample collection, process tests, upload and release results, and manage everything around them — the test catalog, prices, discounts, loyalty rewards, service zones, advertisements and reports.

Who should read this manual

This manual is written for every staff member who signs in to the dashboard. Each role sees only the modules an administrator has enabled for it, so some chapters may describe pages you cannot open.

Role	Typical responsibilities
Admin	Full access to every module. Manages staff accounts, role permissions and system settings.
Manager	Oversees operations, pricing rules, promotions and reports.
Reception	Creates and updates orders, records payments, manages customer accounts.
Lab technician	Processes samples, uploads result PDFs, runs AI review and releases results.
Collector	Field staff who pick up samples. Assigned to collection routes; mainly uses the mobile app.

Note: The exact modules each role can open are configurable. If a page shows *"You don't have permission to access this page"*, ask an administrator to enable it for your role (see chapter 19, Permissions).

How to read this manual

- Button, menu and field names appear in bold, exactly as they appear on screen — for example **Create order**.
- Numbered lists are procedures: follow the steps in order.
- Green **Tip** boxes give shortcuts; blue **Note** boxes add context; red **Important** boxes warn about actions that affect patients or cannot be undone.
- All amounts are in Myanmar Kyat (MMK).

2. Getting Started

Signing in, finding your way around, and personal preferences.

Signing in

1. Open the Lab Admin Web address in a modern browser (Chrome, Edge, Firefox or Safari).
2. On the **Staff sign in** screen, enter your lab staff **Email** and **Password**. Use the eye icon to show or hide the password while typing.
3. Tick **Keep me signed in** if you are on a private, trusted computer.
4. Click **Sign in**. You land on the first module your role can access (normally **Orders**).

Note: Forgot your password? There is no self-service reset — contact your administrator, who can set a new password for you from the Staff page.

The screen layout

After signing in, every page shares the same frame:

Area	What it contains
Sidebar (left)	The IDHC logo and one link per module you are allowed to use, with Logout at the bottom. Click the logo to collapse the sidebar to icons only; click again to expand. Your choice is remembered on this computer.
Header (top)	The current page title, your profile button (photo or initials, name and role), the language menu, the theme menu and the notification bell.
Content area	The page itself: a banner describing the page, filters, the main table and action buttons.

On narrow screens (tablets and phones) the sidebar is hidden. Tap the logo at the top-left to open it as a drawer; tap outside it or press **Esc** to close it.

Your profile and password

1. Click your name in the header to open **My profile**.
2. Update your name or email if needed.
3. To change your password, enter your **Current password**, then a **New password** and **Confirm password**. Leave the password fields empty to keep your current password.
4. Click **Save changes**.

Language and theme

- **Language** — the globe menu switches the dashboard between **English** and **Myanmar (Burmese)**.
- **Theme** — choose **Light** (default bright workspace), **Dark** (reduced glare for low-light rooms) or **IDHC** (red & gold brand).

Both settings are personal and stored on your device. They do not change what patients see in the mobile app (that is configured in System settings).

Notifications

The bell icon shows a badge when new notifications arrive (for example, new orders from the mobile app). Click it to open the inbox.

Working with lists and tables

Most pages follow the same pattern, so once you learn one you know them all:

- **Search and filters** sit above the table. Type in the search box and pick filter values, then click **Clear filters** to reset.
- **Pagination** sits below the table: **Previous / Next** and **Rows per page**. The footer shows e.g. *Showing 1–20 of 134*.
- **Row actions** — the three-dot button at the end of each row opens its actions menu (**Edit**, **Delete** and page-specific actions).
- **Refresh** reloads the data from the server.
- **Forms** open in a pop-up window. Required fields are validated when you save; errors appear in red under the field.
- **Confirmations** — deleting anything always asks you to confirm first.
- **Messages** — a short message (toast) confirms success (green) or explains a failure (red).


Signing out

Click **Logout** at the bottom of the sidebar. Always sign out on shared computers.

3. The Order Lifecycle

How an order moves through the lab, and which page handles each stage.

Every order passes through six statuses. Understanding this flow makes the rest of the system easy to follow:



Status	Meaning	What moves it forward
Pending	New order from the mobile app or created by staff. Tests may still need to be assigned.	Assign a collector route on the Collection page → Scheduled.
Scheduled	A collector and pickup time are assigned.	Start collecting on the Collection page → Collecting.
Collecting	The collector is on the way / taking the sample.	Start lab processing on Lab results → Running.
Running	The sample is in the lab and being tested.	Mark lab complete on Lab results → Completed.
Completed	Testing finished; result PDFs can be uploaded and reviewed.	Upload PDFs, run AI review, then Release / Mark delivered → Delivered.
Delivered	Results are released to the patient (in the app and/or as a hard copy).	Final status.

Priority

Each order is either **Urgent** or **Elective**. Urgent orders are highlighted in lists and counted separately in reports — handle them first.

Report delivery method

Method	What the patient receives
Soft copy	Result PDFs released in the patient mobile app.
Hard copy	A printed report handed over in person. PDFs are kept internal and not shared in the app.
Soft & hard copy	Digital PDFs in the app plus a printed copy.

Payment status

Payment is tracked independently of the order status. An order can be **Unpaid**, **Partial received**, **Received** (paid in full), **Pending** or **Failed**. Payments are recorded on the Orders page.

4. Orders

Review orders from the mobile app: confirm the prescription, assign lab tests, record payments and follow each order through to delivery.

The order list

The table lists every order with patient, priority, status, tests, prices, payment progress and the date it was created. Useful columns:

- **Rx** — shows whether a prescription file was uploaded with the order.
- **Tests** — number of tests assigned. Use the **Tests** filter (*Assigned / Missing*) to find orders that still need tests.
- **Original / Final** — price before and after discounts and fees.
- **Paid / Due** — sum of payments recorded, and the amount still owed (with a % paid indicator).
- **Referral fee** — commission generated by this order's tests.

Search by patient, and filter by status, priority and test assignment.

Creating an order

Use this when a patient orders by phone or at the counter.

1. Click **Create order**.
2. Choose the **Ordering user** (the customer account). Create the user first on the Users page if needed.
3. Enter the **Patient name**, **Age**, **Priority** and an optional **Description** (clinical notes).
4. Set the **Report delivery** method.
5. Set the collection location on the map. The system checks which service zone it falls in and applies that zone's **Service fee**. A location outside every zone cannot be priced.
6. Open **Test**, search by name or code, and tick one or more tests. Click **Done**. Selected tests appear in the table below; use **Delete** on a line to remove it.
7. Optionally choose a **Preferred collector**; otherwise the lab assigns one later.
8. Check the **Order pricing** summary (tests, material fee, service fee, estimated amount due) and click **Create**.

Tip: The catalog can be large. The test picker shows a limited number of matches at a time — type part of the name or code to narrow the list.

Viewing order details

Choose **Detail** from the row's actions menu. The detail window shows:

- Patient, ordering user, contact, address, priority, status and report delivery method.
- The **Prescription** image or PDF, with **Open in new tab** and **Download**.
- Pricing breakdown: original price, discount %, material fee, service fee, referral fee and final price.
- Payments: each payment's amount, method, reference and time, plus any loyalty points redeemed.
- Schedule: collecting person, collection time, running time and report-out time.
- Line items: each test with quantity, unit price and subtotal.

Adding tests to a pending order

Orders from the mobile app often arrive with a prescription but no tests. To assign them:

1. Find the order (tip: filter **Tests = Missing**) and choose **Add test** from the actions menu.
2. Read the clinical notes and prescription shown at the top.
3. Open **Test**, search and tick the tests, then click **Add selected tests**.

Note: Tests can only be added or changed while the order is **Pending**. Use **Update order** to edit patient details, priority, delivery method or the test list of a pending order.

Recording a payment

1. Choose **Add payment** from the row's actions menu. The summary shows the order total, already paid and balance remaining.
2. Enter the **Payment amount** — or click **Pay remaining balance** to fill it in.
3. Choose the **Payment method: Cash, Bank transfer** or **Mobile pay**, and add an optional **Reference no.**
4. If the customer has loyalty points, you can **Redeem points**; the MMK value is deducted from the amount due. Points alone may cover the full balance.
5. Review the *After this payment* preview and click **Record payment**.

Deleting an order

Choose **Delete** from the actions menu and confirm. Only delete orders created by mistake.

Important: Deleting an order cannot be undone.

5. Collection & Routing

Plan efficient sample-pickup routes, assign collectors and start collection.

The Collection page has two views: **pending** orders waiting to be routed, and **scheduled** collections already assigned to a collector.

Important: Routes start from the lab location set in **System settings** → **Location & contact**. Set it before planning your first route. You also need at least one staff member with the **Collector** role.

Planning a route

1. In the pending list, tick one or more orders that need a sample pickup.
2. In **Route setup**, confirm the route start (lab location) and choose the **Start time**.
3. Enter the **Collection duration** — minutes spent at each stop.
4. Choose the **Number of collectors**. The stops will be split across that many routes.
5. Click **Plan route**. The system (AI-assisted) orders the stops efficiently and estimates arrival and finish times for each route.

Assigning collectors

1. For each planned route, choose the collector from the list.
2. Adjust individual stop times (**Arrive, From, Until**) if needed.
3. Click **Assign route** (or **Assign all routes**) and confirm. The orders move from **Pending** to **Scheduled** and the collectors see the stops in their app.

Updating a scheduled pickup

1. Switch to the scheduled view and select one order.
2. Click **Update schedule** to change the collector, **Pickup date** or stop times (**Collect from** / **Collect until**), then **Save schedule**.

Starting collection

When the collector sets off for a pickup, select the scheduled order and click **Start collecting**. The order status changes to **Collecting**.

6. Lab Results

Advance orders through lab processing, upload result PDFs, run AI review, and release results to the patient.

Selecting an order

The table shows orders with their status, priority, report delivery method, collector and collect time. Select one order; its **Result entry** panel opens below the table.

Step 1 — Receive the sample

When the sample arrives at the lab (status **Collecting**), click **Start lab processing**. The order becomes **Running**.

Step 2 — Finish lab processing

When all tests are finished, click **Mark lab complete**. The order becomes **Completed** and the result workflow appears: **Upload PDFs** → **AI review** → **Release**.

Step 3 — Upload result PDFs

- Each test on the order has its own card. Click **Upload PDF** to attach the result file for that test, or **Replace PDF** to swap it.
- When one PDF covers several tests, tick those tests and use **Upload to selected** — they are grouped as a **Combined report**.
- Use **View PDF** / **Download PDF** to check an uploaded file.

Step 4 — Run AI review

Run **AI review** on each uploaded PDF. The analysis usually takes a few seconds and returns **AI Analysis Results**, including a *Profile Alignment* check (does the report match the patient?) and a *Medical Verification* summary. A review flagged with an issue should be checked by a technician; fix the PDF and re-run the review if needed.

Note: Every test needs a PDF and a completed AI review before results can be released. If AI review is not configured, release is blocked — contact your administrator.

Step 5 — Release or deliver

Delivery method	Final action
Soft copy	Click Release / Release to patient . The patient is notified and can view the PDFs in the app.
Hard copy	Print the report, set the Handover time and click Save time (the patient is notified), then click Mark delivered after handing it over.
Soft & hard copy	Release the PDFs in the app and also save the handover time for the printed copy.

The order becomes **Delivered**. You can still replace PDFs or re-run AI review afterwards if a correction is needed.

Important: Releasing sends results to the patient. Double-check the PDFs and AI review before you click Release.

7. Lab Tests (Test Catalog)

Maintain the list of tests the lab offers: names, codes, categories, prices and availability.

The catalog table shows each test's name, code, category, base price, discount, discounted price and status. Search by name or code and filter by **Active only** / **Inactive only**.

Adding a test

1. Click **Create lab test**.
2. Enter the **Test name** and an optional **Description** (panel contents, fasting instructions, etc.).
3. Optionally enter a **Test code** (generated from the name if left empty) and a **Category** such as *Chemistry*.
4. Enter the **Base price** in MMK.
5. Leave the test active and click **Create test**.

Editing, deactivating and deleting

- Choose **Edit** from the actions menu to change any field.
- Switch a test to **inactive** to hide it from the mobile app while keeping it in the catalog and in past orders. This is the safe alternative to deleting.
- **Delete** removes the test permanently.

Tip: Discounts are set on the Discounts page, not on the test itself. The catalog just displays the resulting discounted price.

8. Users

Manage patient and partner accounts: contact details, location and loyalty balance.

The **User directory** lists customers with phone, status, address, map coordinates and loyalty point balance. Customers normally register themselves in the mobile app; use this page to help them or to create an account for a walk-in customer.

Adding a user

1. Click **Add user**.
2. Enter the **Name** and **Phone** in international format (+959...). Each phone number can only be registered once.
3. Enter the **Address** and set the location on the map — search for an address, click the map, or use **Use my location**.
4. Set an **Initial password** the customer will use to sign in to the app.
5. Click **Create user**.

Editing a user

Choose **Edit** from the actions menu. Leave the password empty to keep the current one, or type a new one to reset it for the customer. You can also correct the **Total points** balance here.

9. Staff

Add and maintain the people who work in the lab. Roles control what they can do.

Adding a staff member

1. Click **Add staff**.
2. Enter the **Name** and **Email** (used to sign in; must be unique).
3. Choose the **Role**: Admin, Manager, Reception, Lab technician or Collector.
4. Optionally add a **Profile image** (PNG, JPG, GIF or WebP). It appears in the staff list and on collector route assignments.
5. Set an **Initial password** and click **Create staff**. Share the password with the person securely.

Editing, deactivating and resetting passwords

- Choose **Edit** from the actions menu to change details or role.
- To reset someone's password, type a new one; leave it empty to keep the current password.
- Turn the account **inactive** when someone leaves — inactive staff cannot sign in until reactivated.
- You cannot delete your own account.

Note: Collectors must exist here (role **Collector**, active) before routes can be assigned on the Collection page.

10. Ratings & Feedback

Review patient ratings together with the order they refer to.

Each row shows the patient, who submitted the rating, contact details, the order and its tests, amount, star rating (1–5), the patient's remark and the date.

- Use the **Rating** filter to show only e.g. 1-star reviews that need follow-up.
- Search across patient, user, order, tests, status and remark text.

Tip: Low ratings are a good trigger for a follow-up call. Rating trends also appear in Reports.

11. Discounts

Set a flat discount percentage for each lab test, with an optional date range.

Discounted price = base price × (1 – discount %). Rules apply to active catalog tests only.

Adding a discount

1. Click **Add discount**.
2. Search the **Lab tests** list and tick one or more tests (use **Select all** for everything visible).
3. Enter the **Discount %** (0–100). The preview shows the after-discount price for each test.
4. Optionally set a **Start date** and **End date** for a limited-time promotion.
5. Click **Create rule** — one rule is created per selected test.

Edit or delete a rule from the actions menu. Turn a rule **inactive** to pause it without losing it.

Note: Membership tier discounts stack on top of the test discount (see chapter 13).

12. Referral Fees

Configure the referral commission paid for each lab test and report on it.

Referral fee = base price × referral %. The page has two tabs.

Settings tab

1. Click **Add referral fee**.
2. Search and tick one or more tests.
3. Enter the **Referral fee %**; the preview shows the fee in MMK per test.
4. Click **Create rule**. Edit, deactivate or delete rules from the actions menu.

Report tab

Pick a **Start date** and **End date** to see every order that generated a referral fee in that range, with a summary such as *"42 orders generated 315,000 MMK in referral fees"*.

13. Membership Tiers

Reward loyal customers with a tier discount based on lifetime spend.

Customers are automatically placed in the highest active tier they qualify for by **lifetime spend** (MMK) — for example Normal, Silver, Gold. This is separate from redeemable loyalty points.

A tier's discount is **added** to the test's own discount, capped at 100% combined. Example: a test with a 10% discount ordered by a Gold customer (5%) gets 15% off.

Adding a tier

1. Click **Add tier**.
2. Enter the **Tier name**, **Minimum spend (MMK)** and **Discount percent**.
3. Click **Create tier**.

If you delete or deactivate a tier, its customers fall back to the next tier they qualify for.

14. Service Geofences

Define coverage zones on the map and the service (home-collection) fee for each.

When an order is created, the collection address is matched against these zones and the zone's service fee is added to the order.

Adding a zone

1. Click **Add geofence** and enter a zone name.
2. Draw the zone as a rectangle: click one corner on the map, then the opposite corner. Alternatively type **From / To** place names and click **Find on map**.
3. Drag the box to move it or drag its corners to resize. Click **Clear zone** to start over.
4. Enter the **Service fee (MMK)**.
5. Set **Coverage** to **Inside** or **Outside**. When zones overlap, an Inside zone wins over an Outside zone — use this for a cheaper city zone inside a wider regional zone.
6. Click **Create geofence**.

Testing with the fee calculator

The **Fee calculator** panel lets you click any point on the map (or search, or click a zone's bounding box in the table) and press **Calculate fee** to see which zone applies and the fee. A point outside all zones shows "*outside all service zones*".

Note: Changing or deleting a zone does not change existing orders — only new ones.

15. Advertisements

Manage promotional banners shown in the patient mobile app.

1. Click **Add advertisement**.
2. Enter a **Title** and optional **Description**.
3. Add a **Banner image**: **Choose image** to upload, or paste an image URL. A wide landscape image (about 1200×400 px) works best.
4. Optionally set **Start** and **End** date & time to schedule the ad.
5. Click **Create advertisement**.

The list shows each ad's schedule state:

State	Meaning
Live	Active and within its schedule — visible in the app now.
Scheduled	Active but the start date is in the future.
Expired	The end date has passed.
Off	Turned inactive — never shown until reactivated.

16. Loyalty Points

Control how customers earn and redeem points, and audit point balances.

Tab	Purpose
Earn rules	How many points a customer earns per MMK spent. Example: <i>1,000 MMK per batch → 1 point</i> . Rules can have optional start/end dates for seasonal campaigns.
Points by user	Live point balance of every customer, searchable by name or phone — for support queries.
Redemption rate	The global value of one point in MMK when redeeming against an order.
Redemption records	Every earn, redeem and manual adjustment, newest first, filterable by type.

Adding an earn rule

1. On **Earn rules**, click **Add rule**.
2. Enter a **Rule name**, the **MMK spend per batch** and the **Points earned per batch**.
3. Optionally set start and end dates, then click **Create rule**.

Setting the redemption rate

On **Redemption rate**, enter the **MMK per point** value and click **Save rate**. Until a rate above 0 is set, staff cannot redeem points on payments.

17. System Settings

Settings that affect the patient mobile app and shared lab details.

Section	What you can set
Mobile app	The theme and default language patients see in the mobile app.
Location & contact	Lab address and map location (the start point of every collection route), contact phone and contact email shown to patients.
Material fee	The default fee for lab materials (tubes, kits, etc.) added to new orders.

1. Change the values you need. A bar at the bottom shows **Unsaved changes**.
2. Click **Save settings**. The bar changes to **All changes saved** with the save time.

Reset to defaults restores the mobile theme, language and contact fields to the application defaults (after confirmation). Your own dashboard theme and language are not affected.

18. Reports & Analytics

Dashboard KPIs, operational queues, staff performance and financial summaries.

Choose the **Reporting period** (number of days) at the top. The page then shows:

Section	What it tells you
Snapshot	Orders, delivered revenue, urgent orders, users, AI check results and average rating at a glance.
Discount impact	Original value vs. final revenue, total discount given and effective discount %.
Charts	Daily orders, monthly revenue, tests by category and customer rating breakdown.
Pending results queue	Open orders and how many hours each has been waiting.
Turnaround time (TAT)	Time from sample collection to report out, per order, with the average.
Staff activity	Actions logged per staff member and when they last acted.
Sample collection performance	Collections per collector and average time to collect.
Latest reviews	Recent ratings split into positive (4–5 stars), neutral (3 stars) and negative (1–2 stars).

19. Permissions (Admin only)

Choose which modules each staff role can open.

The **Module access by role** matrix has one row per module and one column per role. Tick a box to give that role access, untick to remove it, then click **Save changes**.

- Hidden modules disappear from the sidebar for that role, and direct links show an access-restricted message.
- **Admin always has full access** and cannot be edited, so an administrator can never be locked out.
- Only admins can see the Permissions page.
- Changes apply the next time affected staff load or refresh the dashboard.

20. Troubleshooting & FAQ

Q: I can't sign in.

Check the email and password (the eye icon helps). If your account is inactive or you forgot the password, ask an administrator.

Q: A menu item I need is missing.

Your role doesn't have access to that module. Ask an admin to enable it on the Permissions page.

Q: "Plan route" says to configure the lab location.

Set the lab address and map pin in System settings → Location & contact, then save.

Q: No collectors to choose from.

Add a staff member with role Collector (and make sure it is active) on the Staff page.

Q: I can't add tests to an order.

Tests can only be changed while the order is Pending.

Q: The service fee can't be calculated when creating an order.

The collection location is outside every active service zone. Move the pin or add a zone.

Q: Release is disabled on Lab results.

Make sure every test has a PDF and a finished AI review. For hard copy, save the handover time first. If AI review isn't configured, contact your admin.

Q: Points can't be redeemed on a payment.

Set a redemption rate above 0 in Loyalty points → Redemption rate.

Q: Data looks out of date.

Click **Refresh** on the page, or reload the browser tab.

Q: A banner says to set VITE_API_BASE_URL.

The website can't reach the server. This is a configuration issue — contact your system administrator or IT support.

Security reminders

- Never share your password; each staff member must have their own account.
- Sign out on shared computers and don't use *Keep me signed in* there.
- Patient data is confidential. All access is logged.